

# Selecting the correct **AeroChamber Plus\*** **Flow-Vu\*** Anti-Static Valved Holding Chamber (VHC) for children & adults

**AeroChamber Plus\* VHC is being replaced with AeroChamber Plus\* Flow-Vu\* Anti-Static VHC from August 2020. Technique should always be taught and assessed when prescribing a new VHC.**



**AeroChamber Plus\* Flow-Vu\* Anti-Static VHC with Medium Mask for Children 1-5 years**



**AeroChamber Plus\* Flow-Vu\* Anti-Static VHC with Youth Mouthpiece 5 years +**



**AeroChamber Plus\* Flow-Vu\* Anti-Static VHC with Small Adult Mask**



**AeroChamber Plus\* Flow-Vu\* Anti-Static VHC with Large Adult Mask**



**AeroChamber Plus\* Flow-Vu\* Anti-Static VHC with Mouthpiece**

**AeroChamber Plus\* Flow-Vu\* Anti-Static VHC should be washed weekly. Please refer to manufacturer's patient information leaflet for washing instructions. The AeroChamber Plus\* Flow-Vu\* Anti-Static VHC should also be replaced annually.**

Further advice on spacers is available from [www.ljf.scot.nhs.uk](http://www.ljf.scot.nhs.uk)

Asthma UK <https://www.asthma.org.uk/advice/inhaler-videos>